




The Centec Group

Centec is a privately owned group of companies. A core business for us is customized systems for the production and distribution of high purity water and water for injection (WFI) for pharmaceutical and biotechnological plants in accordance with GMP and FDA guidelines. We are an experienced partner for pure steam generation, CIP- and SIP-processes as well as for services related to product preparation. Centec technology includes a range of high precision process sensors for accurately measuring critical product properties such as the concentration of acidic and caustic solutions and O_2 content. The largest pharmaceutical groups in the world are among our key customers.



Accuracy. Reliability. Centec.

Germany

Centec GmbH Wilhelm-Röntgen-Strasse 10 63477 Maintal Tel.: +49 6181 18 78 0 Fax: +49 6181 18 78 50 info@centec.de

Czech Republic

Centec automatika s.r.o. Pekařská 8/601 155 00 Praha 5 Tel.: +420 257 084 111 Fax: +420 235 518 701 prodej@centec.cz

USA

Centec LLC
P. O. Box 820
Germantown, WI 53022-0820
Tel.: +1 262 251 8209
Fax: +1 262 251 8376
info@centec-usa.com

UK

Centec UK Stalworths, The Street Great Tey, Colchester, Essex, CO6 1JS Tel.: +44 1206 21 19 21 Fax: +44 1206 21 19 16 info@centec-uk.com

Serbia

Centec Serbia
Bogdana Žerajića 34/III
11000 Beograd
Tel.: + 381 11 358 11 24
Fax: + 381 11 358 11 24
info@centec.rs

India

Centec RRR Systems & Sensors Pvt Ltd RRR House, Plot 80, Sector 23 Turbhe Naka, Navi Mumbai - 400 705 Tel.: +91 22 2783 3655 & 2783 1348 Fax: +91 22 2783 4814 mail@centecrrr.com

Brazi

Centec América Latina Ltda Rua Mexico 148 conj. 1004 Centro 20031 142 Rio de Janeiro Tel.: +55 21 2223 2066 centeclatina@terra.com.br



968

Membrane Deaeration

Centec Pharma Systems



Automated
process skids
and high precision
sensors from a
single source.
Centec.

More information at www.centec.de

DGS

The Principle


The Centec membrane deaeration system DGS is a modular skid for the removal of oxygen from water and low viscosity liquids. Using highly effective membrane contactors, extraordinarily low values of residual O_2 are achieved. The presence of even low oxygen levels can increase microbial activity and adversely affect quality and shelf life of pharmaceutical products. Each contactor contains thousands of microporous hydrophobic hollow fibre membranes. Their large internal surface maximizes the contact area between liquid and gas. While the strip gas $(CO_2 \text{ or } N_2)$ is applied on the inside of the hollow fibres, the water flows in counter current on the outside. The high partial pressure difference of O_2 forces the oxygen out of the liquid to permeate through the membranes and away with the strip gas. This fundamental scientific principle is described by "Henry's Law". The DGS system consumes minimal energy and strip gas. Due to its modular design, DGS can easily be adapted to the required capacity (by parallel installation of membranes) and to the required residual oxygen content (by serial installation). O_2 content monitoring can be achieved with highly accurate Centec OXYTRANS optical sensing technology.

Technical Data

Capacity	1 - 500 m³/h
Residual Oxygen	< 10 ppb
Pressure of Operation	0 - 6 bar
Temperature of Operation	2 - 45 °C
Temperature of CIP	max. 85 °C
Material	1.4404/1.4435/ AISI 316L/ electropolishing possible
Membrane	polypropylene (FDA approved)
PLC	SIMATIC S7
Options	in-line O₂ measurement chilling/chilling water recovery pre-filtration disinfection

The Centec production is certified according to ISO 9001.





- Modular Design with Standard PLC
 skid mounted for easy installation and start-up
- In Compliance with USP, ASME, GMP, FDA, ISPE
- Completely Qualified (DQ, IQ, OQ)
- Low Operating and Maintenance Costs

Experience. Expertise. Centec.

Particle Pre-Filtration · Disinfection · Water Softening & Demineralization · Ultrafiltration · Reverse Osmosis

Electro Deionization · WFI Distillation · Membrane Deaeration · Column Deaeration · Vacuum Deaeration

Multi Component Mixing · Additive Dosing · Flash Pasteurization · Cleaning-in-Place · Pure Steam Generation

