




The Centec Group

Centec is a privately owned group of companies. A core business for us is customized systems for the production and distribution of high purity water and water for injection (WFI) for pharmaceutical and biotechnological plants in accordance with GMP and FDA guidelines. We are an experienced partner for pure steam generation, CIP- and SIP-processes as well as for services related to product preparation. Centec technology includes a range of high precision process sensors for accurately measuring critical product properties such as the concentration of acidic and caustic solutions and O_2 content. The largest pharmaceutical groups in the world are among our key customers.





Germany

Centec GmbH Wilhelm-Röntgen-Strasse 10 63477 Maintal Tel.: +49 6181 18 78 0 Fax: +49 6181 18 78 50 info@centec.de

Czech Republic

Centec automatika s.r.o. Pekařská 8/601 155 00 Praha 5 Tel.: +420 257 084 111 Fax: +420 235 518 701 prodej@centec.cz

USA

Centec LLC
P. O. Box 820
Germantown, WI 53022-0820
Tel.: +1 262 251 8209
Fax: +1 262 251 8376
info@centec-usa.com

UK

Centec UK Stalworths, The Street Great Tey, Colchester, Essex, CO6 1JS Tel.: +44 1206 21 19 21 Fax: +44 1206 21 19 16 info@centec-uk.com

Serbia

Centec Serbia
Bogdana Žerajića 34/III
11000 Beograd
Tel.: + 381 11 358 11 24
Fax: + 381 11 358 11 24
info@centec.rs

India

Centec RRR Systems & Sensors Pvt Ltd RRR House, Plot 80, Sector 23 Turbhe Naka, Navi Mumbai - 400 705 Tel.: +91 22 2783 3655 & 2783 1348 Fax: +91 22 2783 4814 mail@centecrrr.com


Brazil

Centec América Latina Ltda Rua Mexico 148 conj. 1004 Centro 20031 142 Rio de Janeiro Tel.: +55 21 2223 2066 centeclatina@terra.com.br

centec

Product Sterile Filtration

Centec Pharma Systems



Automated
process skids
and high precision
sensors from a
single source.
Centec.

More information at www.centec.de

INTRAP


The Principle

The Centec trap filter skid INTRAP is a single- or multi-stage system containing back flushable filter cartridges with specific pore sizes. This unit can be used for many applications, e.g. for sterile product filtration or as an early step in a water treatment regime. By sterile filtration, particles and micro-organisms are removed from the product. This leads to an improvement of product stability and shelf life. When water is purified, prior solids removal avoids clogging and damage to the membranes that subsequently achieve even higher purity. INTRAP can be designed for manual operation with pressure indicators or as a fully automated system. For automatic back flush regeneration the differential pressure needs to be continuously measured. With this option included, the lifetime of the filter cartridges is maximized. Depending on product characteristics, in some applications sterile filtration can be an alternative to flash pasteurization, which is the elimination of the most harmful micro-organisms by using heat. Other processors may prefer to use sterile filters to treat the product before they flash pasteurize.

Technical Data

Capacity	5 - 150 m ³ /h
Pressure of Operation	2 - 8 bar
Temperature of Operation	2 - 85 °C
Temperature of CIP	up to 85 °
Material	1.4404/1.4435/ AISI 316L/ electropolishing possible
Surface Finish	Ra < 0,4 possible
Ferrite Content	< 1 % possible
PLC	SIMATIC S7
Options	in-line O_2 measurement pressure measurement flash pasteurization





- Modular Design with Standard PLC
 skid mounted for easy installation and start-up
- In Compliance with USP, ASME, GMP, FDA, ISPE
- Completely Qualified (DQ, IQ, OQ)
- Low Operating and Maintenance Costs

Experience. Expertise. Centec.

Particle Pre-Filtration · Disinfection · Water Softening & Demineralization · Ultrafiltration · Reverse Osmosis

Electro Deionization · WFI Distillation · Membrane Deaeration · Column Deaeration · Vacuum Deaeration

Multi Component Mixing · Additive Dosing · Flash Pasteurization · Cleaning-in-Place · Pure Steam Generation

